

**Appendix E: Summary of Safety and Effectiveness Data*****I. General Information***

Company : Fotona d.d.  
Stegne 7, 1210 Ljubljana, Slovenia

Contact Person : Mojca Valjavec

Preparation Date : 01-10-00

Device Trade Name : Depilase YAGLASE Nd:YAG Laser system and Accessories

Common Name : Nd:YAG Pulsed Surgical Laser System

Classification Name : Instrument, Surgical, Powered, Laser  
79-GEX  
21 CFR 878-48

***II. Description***

The Depilase YAGLASE system is based on Nd:YAG laser technology. Within the system, an optical cavity contains the Nd:YAG crystal, which is activated by means of the use of flashlamps. After the cavity, a red diode aiming beam is reflected onto a coaxial beam path using a beamsplitter assembly. The combined therapeutic and aiming beams are guided down an optical fiber delivery system to a focusing handpiece. The laser is used in non-contact mode.

The YAGLASE system is designed with 5 major sub-systems:

- a) A high voltage power supply which converts and rectifies the a.c. mains current to provide regulated power for the flashlamp simmer current and main triggering pulse.
- b) A cooling system consisting of an internal water flow circuit together with water-to-air heat exchanger.
- c) An Nd:YAG laser rod, capable of generating optical pulses at a frequency up to 10 Hz.
- d) An optical delivery system, interfacing the energy from the laser to the patient via an optical fiber and focusing handpiece.
- e) The microprocessor based controller which regulates the functions of the laser and allows parameter selection by the user.

## *II. Intended Use*

The Depilase YAGLASE Nd:YAG Laser System is indicated for the coagulation and hemostasis of vascular lesions and soft tissue

## *III. Summary of Substantial Equivalence*

Depilase believes that its YAGLASE system is substantially equivalent to the Laserscope Lyra (K990903), HGM VeinLase Nd:YAG Laser (K981952), and other Nd:YAG Lasers previously cleared for coagulation and hemostasis of vascular lesions and soft tissue.

All above mentioned lasers are cleared for coagulation and hemostasis of vascular lesions and soft tissue. They therefore have the same Intended Use as the Depilase YAGLASE Laser System..

Technologically, the predicate devices have identical characteristics to the YAGLASE, all comprising a flashlamp pumped Nd:YAG laser rod generating light at a wavelength of 1064 nm, which is subsequently delivered to the patient via an optical fiber delivery system and focusing handpiece.

The YAGLASE Laser output characteristics are very similar to those of predicate devices.

All lasers are microprocessor controlled devices.

All lasers utilize class I aiming beams which pose no hazard to the user.

All systems utilize an internal closed loop water-air heat exchanger circuit for optimal thermal control of the laser cavity.

The risk and benefits for the Depilase YAGLASE are comparable to the predicate devices when used for similar clinical applications.

It is therefore believed that there are no new questions of Safety or Effectiveness raised by the introduction of this device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

APR 18 2000

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Mojca Valjavec, Dipl. Eng.  
Fotona d.d.  
Stegne 7, 1210 Ljubljana  
Slovenia

Re: K000106  
Trade Name: Depilase YAGLASE Nd:YAG Laser System  
and Accessories  
Regulatory Class: II  
Product Code: GEX  
Dated: February 18, 2000  
Received: February 24, 2000

Dear Dr. Valjavec:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Mojca Valjavec, Dipl. Eng.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4595. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



*for* Celia M. Witten, Ph.D., M.D.  
Director  
Division of General and  
Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## INDICATIONS FOR USE STATEMENT

510(k) Number (if known): K000106Device Name: **Depilase YAGLASE Nd:YAG LASER SYSTEM**

## Indications For Use:

The Depilase YAGLASE Nd:YAG Laser System and Accessories are intended for incision, ablation, vaporization, and coagulation and hemostasis of vascular lesions and soft tissue in various surgical areas.

**Dermatology:** Photocoagulation of pigmented lesions to reduce lesion size

**Plastic Surgery:** Coagulation and vaporization of soft tissue.

**Otorhinolaryngology / Head and Neck (ENT):** Tissue ablation and vessel hemostasis.

**Hemostasis during Surgery:** Adjunctive coagulation and hemostasis (bleeding control) during surgery in endoscopic (e.g. laparoscopic) and open procedures.

**Orthopedics:** Ablation, vaporization, incision, excision, coagulation, and hemostasis of soft and cartilaginous tissue in small and large joints including but not limited to knee meniscectomy, knee synovectomy, chondromalacia and tears, shoulder debridement of scar tissue, and synovectomy of the shoulder.

**Neurosurgery:** Hemostasis in neurosurgery procedures such as excision of brain lesions, spinal cord lesions, cranial nerves, peripheral nerves, and pituitary glands.

**Gastroenterology:** Tissue ablation and hemostasis in the gastrointestinal tract; esophageal neoplastic obstructions including squamous cell carcinoma and adenocarcinoma;  
Gastrointestinal hemostasis including varices, esophagitis, esophageal ulcer, Mallory-Weiss tear, gastric ulcer, stomach ulcers, angiodysplasia, non-bleeding ulcers, gastric erosions;  
Gastrointestinal tissue ablation including benign and malignant neoplasm, angiodysplasia, polyps, ulcer, colitis, hemorrhoids.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Diana R. Lochner.  
(Division Sign-Off)  
Division of General Restorative Devices  
510(k) Number K000106

Prescription Use X  
(Per 21 CFR 801.109)

OR

Over-The-Counter Use \_\_\_\_\_